# Séquençage

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Séquence TEM

# http://web.expasy.org/translate/

The screenshot shows a Windows Internet Explorer browser window with the title "ExPASy - Translate tool - Windows Internet Explorer". The address bar contains the URL "http://web.expasy.org/translate/". The browser's menu bar includes "Fichier", "Edition", "Affichage", "Favoris", and "Outils". The search bar features the Google logo and the text "Rechercher" and "Autres >>". The favorites bar shows "Favoris" and several links, including "Aurore boréale photos, vid...", "Sites suggérés", "CHU de Liège Bienvenue", "Accéder à plus de mo...", and "Hotmail". The main content area has a header with the ExPASy logo (SIB and ExPASy Bioinformatics Resource Portal) and the word "Translate". Below the header, the text reads "Translate tool" followed by a description: "Translate is a tool which allows the translation of a nucleotide (DNA/RNA) sequence to a protein sequence." A prompt asks the user to "Please enter a DNA or RNA sequence in the box below (numbers and blanks are ignored)." Below this is a large, empty text input box with a vertical scrollbar. At the bottom, there is a dropdown menu for "Output format:" set to "Verbose ('Met', 'Stop', spaces between residues)", and two buttons: "Reset" and "TRANSLATE SEQUENCE".

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SIB ExpASY Bioinformatics Resource Portal Translate

**Translate tool**

Translate is a tool which allows the translation of a nucleotide (DNA/RNA) sequence to a protein sequence.

Please enter a DNA or RNA sequence in the box below (numbers and blanks are ignored).

```
CATGATAATAATGGTTTCTTAGACGTCAGGTGGCACTTTTCGGGGAAATGTGCGCGGAACCCCTATTTGT
TTATTTTTCTAAATACATTCAAATATGATCCGCTCATGAGACAATAACCCCTGATAAATGCTTCAATAAT
ATTGAAAAAGGAAGAGTATGAGTATTC AACATTTTCGTGTGCGCCCTTATCCCTTTTTTGCGGCATTTTG
CCTTCCTGTTTTGCTCACCCAGAAACGCTGGTGAAGTAAAAGATGCTGAAGATCAGTTGGGTGCACGA
GTGGGTACATCGAACTGGATCTCAACAGCGGTAAGATCCTTGAGAGTTTTCGCCCCGAAGAACGTTTTTC
CAATGATGAGCACTTTTAAAGTTCTGCTATGTGGTGGGTATTATCCCGTGTGACGCCGGGCAAGAGCA
ACTCGGTGCGCGCATACACTATTCTCAGAATGACTTGGTTGAGTACTCACCAAGTACAGAAAAGCATCTT
ACGGATGGCATGACAGTAAGAGAATTATGCAATGCTGCCATACCATGAGTGATAACACTGCTGCCAAT
TACTTCTGACAACGATCGGAGGACCGAAGGAGCTAACCGCTTTTTTGACACAACATGGGGGATCATGTAAC
TCGCCCTTGATCGTTGGGAACCGGAGCTGAATGAAGCCATACCAAACGACGAGCGTGACACCACGATGCCT
GCAGCAATGGCAACAACGTTGCGCAAACTATTAACGGGCAACTACTTACTCTAGCTTCCCGGCAACAAT
TAATAGACTGGATGGAGGCGGATAAAGTTGCAGGACCCTTCTGCGCTCGGCCCTCCGGCTGGCTGGTT
TATTGCTGATAAATCTGGAGCCGGTGAGCGTGGGTCTCGCGGTATCATTGCAGCACTGGGGCCAGATGGT
AAGCCCTCCGCTATCGTAGTTATCTACACGACGGGGAGTCAGGCAACTATGGATGAACGAAATAGACAGA
TCGCTG
```

Output format:

or

Please select one of the following frames:

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http://web.expasy.org/cgi-bin/translate/dna\_sequences?/work/expasy/tmp/http/seqdna.7701,2

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SIB ExpASY Bioinformatics Resource Portal Translate Home Contact

### Translate Tool

Please select one of the "Methionine" or one of the highlighted residues following a **Stop** codon (or the beginning of the sequence).

This will create a virtual Swiss-Prot entry, comprising the residues from your chosen start position up to the following **Stop** codon.

MIIMVS **Stop** TSGGTFRGNVRGTPICLFF **Stop** IHSNMYPLMRQ **Stop** P **Stop** \_ MLQ **Stop** Y **Stop** KRK**SMS** QHFRVALIPFFAAFCLPVFAH  
PETLVKVKDAEDQLGARVGYIELDLNSGKILESFRPEERFPMMSTFKVLLCGAVLSRVDAQQEQLGRRIHYSQNDLVEYSPVTEK  
HLTDG **M**TVRELCSAAIT **M**SDNTAANLLLTIGGPKELTAFLHN **M**GDHVTRLDRWEPNELNEAIPNDERDTT **M**PAAMATTLRKLTTGE  
LLTLASRQQLIDW **M**EADKVVAGPLLRSALPAGWFIADKSGAGERGSRGIIAALGPDGKPSRIVVIYTTGSQAT **M**DERNRQIA

# http://web.expasy.org/translate/

ExPASy Translate Tool - Virtual entry - Windows Internet Explorer

http://web.expasy.org/cgi-bin/translate/dna\_sequences?/work/expasy/tmp/http/seqdna.7701,2,53

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Sequence in FASTA format

[BLAST](#) BLAST submission on ExPASy/SIB

 Sequence analysis tools: ProtParam, ProtScale, Compute pI/Mw, PeptideMass, PeptideCutter,

 ScanProsite

 Direct Submission to SWISS-MODEL

# http://web.expasy.org/translate/

ExPASy Translate Tool - Virtual entry - Windows Internet Explorer

http://web.expasy.org/cgi-bin/translate/dna\_sequences?work/expasy/tmp/http/seqdna.7701,2,53

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Direct Submission to SWISS-MODEL

# http://www.ncbi.nlm.nih.gov/blast/

The screenshot shows the NCBI BLAST web interface in a Windows Internet Explorer browser. The browser's address bar displays the URL: [http://www.ncbi.nlm.nih.gov/blast/Blast.cgi?PAGE=Proteins&PROGRAM=blastp&BLAST\\_PROGRAMS=blastp&PAGE\\_TYPE=BlastSearch&SHOW\\_DEFAULTS=on](http://www.ncbi.nlm.nih.gov/blast/Blast.cgi?PAGE=Proteins&PROGRAM=blastp&BLAST_PROGRAMS=blastp&PAGE_TYPE=BlastSearch&SHOW_DEFAULTS=on). The browser's menu bar includes 'Fichier', 'Edition', 'Affichage', 'Favoris', and 'Outils'. The search bar contains the text 'Rechercher' and 'Autres >>'. The browser's toolbar shows several open tabs, including 'ExpASY Translate Tool - Virtu...' and 'Protein BLAST: search pr...'. The BLAST interface is titled 'Standard Protein BLAST' and features a navigation bar with 'Home', 'Recent Results', 'Saved Strategies', and 'Help'. The 'blastp' tab is highlighted with a red circle. The main content area includes a 'Query subrange' section with 'From' and 'To' input fields, and a 'Choose Search Set' section with a 'Database' dropdown set to 'Non-redundant protein sequences (nr)'. The 'Job Title' field is empty, and the 'Align two or more sequences' checkbox is unchecked.

# http://www.ncbi.nlm.nih.gov/blast/

Protein BLAST: search protein databases using a protein query - Windows Internet Explorer

http://www.ncbi.nlm.nih.gov/blast/Blast.cgi?PAGE=Proteins&PROGRAM=blastp&BLAST\_PROGRAMS=blastp&PAGE\_TYPE=BlastSearch&SHOW\_DEFAULTS=on

Home Recent Results Saved Strategies Help

NCBI/ BLAST/ blastp suite **Standard Protein BLAST**

blastn **blastp** blastx tblastn tblastx

BLASTP programs search protein databases using a protein query. [more...](#) [Reset page](#) [Bookmark](#)

Enter Query Sequence

Enter accession number(s), gi(s), or FASTA sequence(s) [Clear](#) Query subrange [From](#) [To](#)

IPNDERDTIM  
PAAAMATTLRK LLIGELLTLA SRQQIDWME ADKVAGPLLR SALPAGWFIA  
DKSGAGERGS  
RGIIAALGPD GKPSRIVVIY TTGSQATMDE RNRQIA

Or, upload file  [Parcourir...](#)

Job Title   
Enter a descriptive title for your BLAST search

Align two or more sequences

PHI-BLAST (Pattern Hit Initiated BLAST)

DELTA-BLAST (Domain Enhanced Lookup Time Accelerated BLAST)

Choose a BLAST algorithm

**BLAST** Search database Non-redundant protein sequences (nr) using Blastp (protein-protein BLAST)

show results in a new window

[Algorithm parameters](#)

BLAST is a registered trademark of the National Library of Medicine.

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**BLAST** Basic Local Alignment Search Tool My NCBI [Sign In] [Register]

Home Recent Results Saved Strategies Help

NCBI/ BLAST/ blastp suite/ Formatting Results - PVD06TAX012 [\[Formatting options\]](#)

**Job Title: Protein Sequence (276 letters)**

Putative conserved domains have been detected, click on the image below for detailed results.

Query seq. 50 100 150 200 250 276

Superfamilies Transpeptidase superfamily

Multi-domains PRK15442

Request ID	PVD06TAX012
Status	Searching
Submitted at	Sat Mar 24 07:46:14 2012
Current time	Sat Mar 24 07:46:16 2012
Time since submission	00:00:02

This page will be automatically updated in 12 seconds

BLAST is a registered trademark of the National Library of Medicine.

[Copyright](#) | [Disclaimer](#) | [Privacy](#) | [Accessibility](#) | [Contact](#) | [Send feedback](#) [NCBI](#) | [NLM](#) | [NIH](#) | [DHHS](#)

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
ExPASy Translate Tool - Virtu... NCBI Blast: Protein Sequ...

**Descriptions**

Legend for links to other resources: UniGene GEO Gene Structure Map Viewer PubChem BioAssay

Sequences producing significant alignments:

Accession	Description	Max score	Total score	Query coverage	E value	Max ident	Links
<a href="#">AAS4684.1</a>	TEM-1 beta-lactamase [Serratia marcescens] >qb AAS46846.1  TEM-1	571	571	100%	0.0	100%	
<a href="#">AAS46845.1</a>	TEM-1 beta-lactamase [Serratia marcescens]	570	570	100%	0.0	100%	
<a href="#">ACT97477.1</a>	TEM-1b beta-lactamase [mixed culture bacterium CA_qf1DD01_19]	571	571	100%	0.0	100%	
<a href="#">ABY81258.1</a>	beta-lactamase [Escherichia coli] >qb ACK75718.1  beta-lactamase	570	570	100%	0.0	100%	
<a href="#">CAL26907.1</a>	beta-lactamase protein [Acetobacter pasteurianus]	570	570	100%	0.0	100%	
<a href="#">ACZ02437.2</a>	extended spectrum beta-lactamase TEM-1 [Providencia rettgeri]	570	570	100%	0.0	100%	
<a href="#">AAQ57123.1</a>	beta-lactamase TEM-128 [Acinetobacter baumannii]	570	570	100%	0.0	100%	
<a href="#">NP_052173.1</a>	beta-lactamase [Neisseria gonorrhoeae] >ref NP_569411.1  beta-lac	570	570	100%	0.0	100%	
<a href="#">ABP57022.1</a>	beta lactamase [Cloning vector pTrueBlue]	569	569	100%	0.0	99%	
<a href="#">AAD13350.1</a>	beta-lactamase [Expression vector pPK113] >qb AAD29959.1  beta-	569	569	100%	0.0	99%	
<a href="#">ACF32746.1</a>	putative inhibitor-resistant TEM-163 beta-lactamase [Escherichia co	569	569	100%	0.0	99%	
<a href="#">AEF15432.1</a>	beta-lactamase TEM-1 [Escherichia coli]	568	568	99%	0.0	100%	
<a href="#">ACN82383.1</a>	TEM-116 beta-lactamase [Aeromonas hydrophila] >qb ACN82384.1	569	569	100%	0.0	99%	
<a href="#">ADL40245.1</a>	TEM-1 beta-lactamase [Pseudomonas aeruginosa] >qb ADR29904.1	569	569	100%	0.0	99%	
<a href="#">AAF74292.1</a>	extended-spectrum beta-lactamase [Escherichia coli]	568	568	100%	0.0	99%	
<a href="#">CAG30723.1</a>	hypothetical protein [Hordeum vulgare subsp. vulgare] >qb ACT9842	568	568	100%	0.0	99%	
<a href="#">YP_003829050.1</a>	beta-lactamase class A [Escherichia coli] >qb ADL13944.1  TEM-33	568	568	100%	0.0	99%	
<a href="#">CAA74912.2</a>	beta-lactamase class A [Capnocytophaga ochracea]	568	568	100%	0.0	99%	
<a href="#">AAA72841.1</a>	b-lactamase [unidentified cloning vector] >qb AFB73755.1  beta-lac	568	568	100%	0.0	99%	
<a href="#">AAM61953.1</a>	beta-lactamase TEM-105 [Escherichia coli]	568	568	100%	0.0	99%	
<a href="#">ACX83572.1</a>	extended spectrum beta-lactamase TEM-1 [Morqanella morqanii]	568	568	99%	0.0	100%	
<a href="#">CAJ66089.1</a>	beta-lactamase TEM-150 [Escherichia coli] >qb ACT32331.1  extend	568	568	100%	0.0	99%	

# 2<sup>ème</sup> exemple

Séquence SHV

# http://web.expasy.org/translate/

The screenshot shows a web browser window titled "ExPASy - Translate tool - Windows Internet Explorer". The address bar displays "http://web.expasy.org/translate/". The browser's menu bar includes "Fichier", "Edition", "Affichage", "Favoris", and "Outils". A search bar with the Google logo is present, along with a search button labeled "Rechercher" and a link "Autres >>". The browser's favorites bar shows several items, including "Aurore boréale photos, vid...", "Sites suggérés", "CHU de Liège Bienvenue", "Accéder à plus de mo...", and "Hotmail". The main content area of the browser displays the ExPASy website interface. At the top, there is a logo for "SIB" and "ExPASy Bioinformatics Resource Portal" on the left, and the word "Translate" on the right. Below this, the heading "Translate tool" is followed by a description: "Translate is a tool which allows the translation of a nucleotide (DNA/RNA) sequence to a protein sequence." A text input box contains a long DNA sequence: "ATGCGTTATATTCGCTGTGTATTATCTCCCTGTTAGCCACCCTGCGGCTGGCGGTACAGCCAGCCCGCAGCCGCTTGAGCAAATTAACAAGCGAAGCCAGCTGTTCGGCCCGCGTAGGCATGATAGAAATGGATCTGGCCAGCGGCCGACGCTGACCGCCTGGCCGCGCGATGAACGCTTCCCATGATGAGCACCTTTAAAGTATGTCTCTGCGGCGCAGTGCTGGCGCGGGTGGATGCCGGTGACGAACAGCTGGAGCGAAAAGATCCACTATCGCCAGCAGGATCTGGTGGACTACTCGCCGGTCAGCGAAAAACACCTTGCCGACGGCATGACGGTTCGGCGA ACTCTGCGCCCGCCGATTACCATGAGCGATAACAGCGCCGCAATCTGCTGCTGGCCACCGTCGGCGGC CCGCAGGATTGACTGCCTTTTGGCCAGATCGGGGACAAACGTCACCCGCTTGACCGCTGGGAAACGG AACTGAATGAGGCGCTTCCGGGGACGCGCGGACACCCTACCCGGCCAGCATGGCCGCGACCCTGGG CAAGCTGCTGACCAGCCAGCGTCTGAGCGCCCGTTCCGCAACGGCAGCTGCTGCACTGGATGGTGGACGAT CGGGTCGCGGACCGTTGATCCGCTCCGTGCTGCCGGCGGGCTGGTTTATCGCCGATAAGACCGGAGCTA GCAAGCGGGGTGCGCGCGGGATTGTCCGCCTGCTTGCCCGGAATAACAAAGCAGAGCGCATTGTGGTGAT TTATCTGCGGGATACGCGCGGAGCATGGCCGAGCGAAATCAGCAAATCGCCGGGATCGGCGCGCGTGA TCGAG". Below the input box, there is a dropdown menu for "Output format" set to "Verbose ('Met', 'Stop', spaces between residues)". At the bottom, there are two buttons: "Reset" and "TRANSLATE SEQUENCE".

# http://web.expasy.org/translate/

Translate Tool - Results of translation - Windows Internet Explorer

http://web.expasy.org/cgi-bin/translate/dna\_aa

Translate

Translate Tool - Results of translation

Please select one of the following frames:

**5'3' Frame 1**  
Met RYIRLCIISV L AT L P L A V H A S P Q P L E Q I K Q S E S Q L S G R V G Met IE Met D L A S G R T L T A W R A D E R F P Met Met S T F K V V L C G A V L A R V  
S A G D E O L E R K I H Y R Q Q D L V D Y S P V S E K H L A D G Met T V G E L C A A A I T Met S D N S A A N L L L A T V G G P A G L T A F L R Q I G D N V T R L D R W E  
T E L N E A L P G D A R D T T T P A S Met A A T L R K L L T S Q R L S A R S Q R Q L L Q W Met V D D R V A G P L I R S V L P A G W F I A D K T G A S K R G A R G I V A L  
L G P N N K A E R I V V I Y L R D T P A S Met A E R N Q Q I A G I G A A S Q R S

**5'3' Frame 2**  
C V I F A C V L S P C Stop P P C R W R Y T P A R S R L S K L N K A K A S C R A A Stop A Stop Stop K W I W P A A A R Stop P P G A P Met N A F P Stop Stop A P L K  
Stop C S A A Q C W R G W Met P V T N S W S E R S T I A S R I W W T T R R S A K N T L P T A Stop R S A N S A P P L P Stop A I T A P P I C C W P P S A A P Q D Stop  
L P F C A R S A T T S P A L T A G K R N Stop Met R R R F P A T P A T P L P R P A W P R P C A S C Stop P A S V Stop A P V R N G S C C S G W W T I G S P D R Stop S A  
P C C R R A G L S P I R P E L A S G V R A G L S P C L A R I T K Q S A L W Stop F I C G I R R R A W P S E I S K S P G S A R R D R

**5'3' Frame 3**  
A L Y S P V Y Y L P V S H P A A G G T R Q P A A A Stop A N Stop T K R K P A V G P R R H D R N G S G Q R P H A D R L A R R Stop T L S H D E H L Stop S S A L R R S A  
G A G G C R Stop R T A G A K D P L S P A G S G G L L A G Q R K T P C R R H D G R R T L R R R H Y H E R Stop Q R R Q S A A G H R R R P R R I D C L F A P D R R Q R  
H P P Stop P L G N G T E Stop G A S R R R P R H H Y P G Q H G R D P A Q A A D Q P A S E R P F A T A A A A V D G G R S G R R T V D P L R A A G G L V Y R R Stop D  
R S Stop Q A G C A R D C R P A W P E Stop Q S R A H C G D L S A G Y A G E H G R A K S A N R R D R R G V I E

**3'5' Frame 1**  
L D H A A P I P A I C Stop F R S A Met L A G V S R R Stop I T T Met R S A L L F G P S R A T I P R A P R L L A P V L S A I N Q P A G S T E R I N G P A T R S S T I H C S S C  
R C E R A L R R W L V S S L R R V A A Met L A G V V V S R A S P G S A F S S V S Q R S R R V T L S P I W R K K A V N P A G P P T V A S S R L A A L L S L Met V Met  
A A A Q S S P T V Met P S A R C F S L T G E Stop S T R S C W R Stop W I F R S S C S S P A S T R A S T A P Q S T T L K V L I Met G K R S S A R Q A V S V R P L A R S I  
S I Met P T R P D S W L S L C I C S S G C L A C T A S G R V A N R E I I H R R I Stop R

**3'5' Frame 2**  
S I T P R R S R R F A D F A R P C S P A Y P A D K S P Q C A L L C Y S G Q A G R Q S R A H P A C Stop L R S Y R R Stop T S P P A A R S G S T V R R P D R P P S T A A  
A A V A N G R S D A G W S A A C A G S R P C W P G Stop W C R G R R R E A P H S V P F P S G Q G G Stop R C R R S G A K R Q S I L R G R R R W P A A D W R R C Y  
R S W Stop W R R R R V R R P S C R R Q G V F R Stop P A S S P D P A G D S G S F A P A V R H R H P P A P A L R R R A L L Stop R C S S W E S V H R R A R R S A  
C G R W P D P F L S C L R G P T A G F R F V Stop F A Q A A A G W R V P P A A G W L T G R Stop Y T G E Y N A

**3'5' Frame 3**  
R S R R A D P G D L L I S L G H A R R R I P Q I N H H N A L C F V I R A K Q G D N P A R T P L A S S G L I G D K P A R R Q H G A D Q R S G D P I V H H P L Q Q L P L R T G  
A Q T L A G Q Q L A Q G R G H A G R G S G V A G K R L I Q F R P A V K A G D V V A D L A Q K G S Q S C G A A D G G Q Q I G G A V I A H G N G G G A E F A D  
R H A V G K V F F A D R R V V H Q I L L A I V D L S L Q L F V T G I H P R Q H C A A E H Y F K G A H H G K A F I G A P G G Q R A A A G Q I H F Y H A Y A A R Q L A F A L F  
N L L K R L R A G V Y R Q R Q G G Stop Q G D N T Q A N I T H

# http://web.expasy.org/translate/

ExPASy Translate Tool - Select an initiator methionine - Windows Internet Explorer

http://web.expasy.org/cgi-bin/translate/dna\_sequences?/work/expasy/tmp/http/seqdna.10965,1

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STB EXPASY Biinformatics Resource Portal Translate Home | Contact

**Translate Tool**

Please select one of the "Methionine" or one of the highlighted residues following a **Stop** codon (or the beginning of the sequence).

This will create a virtual Swiss-Prot entry, comprising the residues from your chosen start position up to the following **Stop** codon.

MRVIRLCIISLLATLPLAVHASPQPLEQIKQSESQLSGRVGMIEMDLASGRTLTAWRADERFPM MSTFKVVLCGAVLARVDAGDE  
QLERKIHRYRQQDLVDYSPVSEKHLADGMTV GELCAAAITMSDNSAANLLLATVGGPAGLTAFLRQIGDNVTRLDRWETELNEALP  
GDARDTTTPASMAATLRKLLTSQRLSARSQRQLLQW MVDDR VAGPLIRS VLPAGWFIADKTGASKR GARGIVALLGPNNKAERIV  
VIYLRDTPASMAERNQQIAGIGAA Stop S

# http://web.expasy.org/translate/

ExPASy Translate Tool - Virtual entry - Windows Internet Explorer

http://web.expasy.org/cgi-bin/translate/dna\_sequences?/work/expasy/tmp/http/seqdna.10965,1,1

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Sequence in FASTA format

[BLAST](#) BLAST submission on ExpASY/SIB

 Sequence analysis tools: ProtParam, ProtScale, Compute pI/Mw, PeptideMass, PeptideCutter,

 ScanProsite

 Direct Submission to SWISS-MODEL

# http://www.ncbi.nlm.nih.gov/blast/

The screenshot shows the NCBI Blast web interface in a Windows Internet Explorer browser. The browser's address bar displays the URL <http://www.ncbi.nlm.nih.gov/blast/Blast.cgi>. The page title is "NCBI Blast:Protein Sequence (279 letters)". The browser's menu bar includes "Fichier", "Edition", "Affichage", "Favoris", and "Outils". The search bar contains the text "Rechercher" and "Autres >>". The browser's toolbar shows various icons for navigation and printing. The main content area displays a table of search results. A red box highlights the first row of the table, which is circled in red. The table has the following columns: Accession, Description, Max score, Total score, Query coverage, E value, Max ident, and Links. The first row is: CAI76\_27.1, beta-lactamase SHV-12 [Escherichia coli], 566, 566, 100%, 0.0, 100%. The second row is: ABD93867.1, SHV-12 extended-spectrum beta-lactamase [Klebsiella pneumoniae], 566, 566, 100%, 0.0, 100%. The third row is: AAS09882.1, SHV-5 extended spectrum beta-lactamase [Escherichia coli], 566, 566, 100%, 0.0, 100%. The fourth row is: YP\_001338660.1, putative beta-lactamase SHV-12 [Klebsiella pneumoniae subsp. pneu], 565, 565, 100%, 0.0, 100%. The fifth row is: AAU84950.1, SHV beta-lactamase [Klebsiella pneumoniae], 565, 565, 100%, 0.0, 99%. The sixth row is: AAR97951.1, beta-lactamase SHV-5 [Escherichia coli], 564, 564, 100%, 0.0, 99%. The seventh row is: ABA06588.1, extended-spectrum beta-lactamase SHV-66 [Klebsiella pneumoniae], 564, 564, 100%, 0.0, 99%. The eighth row is: ADM25824.1, extended spectrum beta-lactamase SHV-134 [Klebsiella pneumoniae], 564, 564, 100%, 0.0, 99%. The ninth row is: AAS09880.1, SHV-5 extended spectrum beta-lactamase [Escherichia coli], 564, 564, 100%, 0.0, 99%. The tenth row is: ABA06586.1, mutant extended-spectrum beta-lactamase SHV-64 [Klebsiella pneu], 564, 564, 100%, 0.0, 99%. The eleventh row is: ACN72587.1, beta-lactamase SHV-12-like protein [Salmonella enterica subsp. entr], 564, 564, 100%, 0.0, 99%. The twelfth row is: POA928.1, RecName: Full=Beta-lactamase SHV-2; AltName: Full=SHV-2A; Flags, 564, 564, 100%, 0.0, 99%. The thirteenth row is: AAS09884.1, SHV-5 extended spectrum beta-lactamase [Escherichia coli], 563, 563, 100%, 0.0, 99%. The fourteenth row is: AAP20890.1, extended-spectrum beta-lactamase SHV-48 [Acinetobacter baumani], 563, 563, 100%, 0.0, 99%. The fifteenth row is: ABC25482.1, beta-lactamase SHV-71 [Acinetobacter baumannii], 563, 563, 100%, 0.0, 99%. The sixteenth row is: AEO09376.1, beta-lactamase, partial [Enterobacter cloacae], 563, 563, 99%, 0.0, 100%. The seventeenth row is: AAP33454.2, extended spectrum beta-lactamase [Klebsiella pneumoniae] >qb|AC, 563, 563, 100%, 0.0, 99%. The eighteenth row is: AAQ55480.1, beta-lactamase SHV-56 [Acinetobacter baumannii], 563, 563, 100%, 0.0, 99%. The nineteenth row is: CAB37325.2, Beta-lactamase [Escherichia coli], 563, 563, 100%, 0.0, 99%. The twentieth row is: ABC58727.1, extended-spectrum beta-lactamase SHV-86 [Klebsiella pneumoniae], 563, 563, 100%, 0.0, 99%.

Legend for links to other resources: [U](#) UniGene [E](#) GEO [G](#) Gene [S](#) Structure [M](#) Map Viewer [B](#) PubChem BioAssay

Sequences producing significant alignments:

Accession	Description	Max score	Total score	Query coverage	E value	Max ident	Links
CAI76_27.1	beta-lactamase SHV-12 [Escherichia coli]	566	566	100%	0.0	100%	
ABD93867.1	SHV-12 extended-spectrum beta-lactamase [Klebsiella pneumoniae]	566	566	100%	0.0	100%	
AAS09882.1	SHV-5 extended spectrum beta-lactamase [Escherichia coli]	566	566	100%	0.0	100%	
YP_001338660.1	putative beta-lactamase SHV-12 [Klebsiella pneumoniae subsp. pneu]	565	565	100%	0.0	100%	<a href="#">G</a>
AAU84950.1	SHV beta-lactamase [Klebsiella pneumoniae]	565	565	100%	0.0	99%	
AAR97951.1	beta-lactamase SHV-5 [Escherichia coli]	564	564	100%	0.0	99%	
ABA06588.1	extended-spectrum beta-lactamase SHV-66 [Klebsiella pneumoniae]	564	564	100%	0.0	99%	
ADM25824.1	extended spectrum beta-lactamase SHV-134 [Klebsiella pneumoniae]	564	564	100%	0.0	99%	
AAS09880.1	SHV-5 extended spectrum beta-lactamase [Escherichia coli]	564	564	100%	0.0	99%	
ABA06586.1	mutant extended-spectrum beta-lactamase SHV-64 [Klebsiella pneu]	564	564	100%	0.0	99%	
ACN72587.1	beta-lactamase SHV-12-like protein [Salmonella enterica subsp. entr]	564	564	100%	0.0	99%	
POA928.1	RecName: Full=Beta-lactamase SHV-2; AltName: Full=SHV-2A; Flags	564	564	100%	0.0	99%	
AAS09884.1	SHV-5 extended spectrum beta-lactamase [Escherichia coli]	563	563	100%	0.0	99%	
AAP20890.1	extended-spectrum beta-lactamase SHV-48 [Acinetobacter baumani]	563	563	100%	0.0	99%	
ABC25482.1	beta-lactamase SHV-71 [Acinetobacter baumannii]	563	563	100%	0.0	99%	
AEO09376.1	beta-lactamase, partial [Enterobacter cloacae]	563	563	99%	0.0	100%	
AAP33454.2	extended spectrum beta-lactamase [Klebsiella pneumoniae] >qb AC	563	563	100%	0.0	99%	
AAQ55480.1	beta-lactamase SHV-56 [Acinetobacter baumannii]	563	563	100%	0.0	99%	
CAB37325.2	Beta-lactamase [Escherichia coli]	563	563	100%	0.0	99%	
ABC58727.1	extended-spectrum beta-lactamase SHV-86 [Klebsiella pneumoniae]	563	563	100%	0.0	99%	

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