

# A simple antibubble maker

S. Dorbolo

*CESAM-GRASP, Département de Physique B5, Université de Liège, B-4000 Liège, Belgium.*

(Dated: November 15, 2019)

PACS numbers:

The system consists in releasing a droplet on a soap film. The height of drop must be sufficient for the droplet to cross the interface. After the passage through the interface, the droplet is wrapped in a bubble. When the wrapped droplet hits the surface of the liquid, the bubble helps the formation of the antibubble. See L. Bai et al, COLSUA **509**, 334 (2016) for more details.

You can use nearly any dishwashing liquid. It is important to work on a liquid surface free of any bubble. Here is the procedure:

1. Prepare soapy water
2. Over fill a tank with the soapy water
3. Fill the syringe with the soapy water and fix the syringe to the antibubble generator
4. Make a soapy film in the frame of the antibubble generator
5. Set the frame about 10-15 mm above the surface of the soapy water

6. Release a droplet.




FIG. 1: (left) The generator is mounted at the tip of a syringe. (right) The generator releasing droplet through the frame. The frame is about 10-15 mm above the surface of the tank. Antibubbles can be seen under the surface.