

# Bidsme – multimodal dataset bidsifier

## Software demonstration

Nikita BELIY<sup>1</sup>

<sup>1</sup>Liege University, Belgium

OSR OHBM, 23 June, 2020



# Introduction

Maintainer: Nikita BELIY [Nikita.Beliy@uliege.be]

## Bidsme

[CRC GitLab](#), [GitHub](#)

- Follows BIDS [v1.2.0](#)
- Allow deviations from BIDS
- MRI (DICOM, Nifti), PET (DICOM, Nifti, ECAT), EEG (BrainVision)
- Extendable  
Data formats and modalities easy to add
- Supports plugins

# Source dataset

[GitHub:datalad/example-dicom-structural](https://github.com/datalad/example-dicom-structural)

- minimal subset from [10.1038/sdata.2014.3 (2014)]
- anonymized, one subject, structural T1-weighted scan
- Modified: for half of dcm files
  - **Series Number** changed to *402*
  - **Series Description** and **Protocol Name** changed to *anat-T2w*

# Step 1: Dataset preparation

```
bidsme prepare dicoms/ prepared/ --no-subject --no-session
```

- **bidsme** – main programm, alias for *python3 bidsme.py*
- **prepare** – command, what the programm will do
- **dicoms/** – source folder, where original dicoms are found
- **prepared/** – destination folder for prepared dataset
- **--no-subject** – dicoms are not stored in folder per-subject
- **--no-session** – dicoms are not stored in folder per-session

## Result:

```
~/tests/bidsme/prepared      File: sub-02
../                          4 K
code/                        4 K
sub-02/                      4 K
participants.json            124 B
participants.tsv              22 B

sub-02/
-- ses-/
-- MRI/
-- 401-anat-T1w/
-- N2D_0001.dcm
-- N2D_0002.dcm
```

## Step 2: Convert to Nifti

```
dcm2niix -b n -f N2D_0001  
prepared/sub-02/ses-/MRI/401-anat-T1w/N2D_0001.dcm
```

- dcm2niix – converter
- **-b n** – do not export meta to json
- **-f N2D\_0001** – conserve the name of file

```
rm prepared/sub-02/ses-/MRI/*/*.dcm
```

- Remove converted dcm files

### Result:

```
prepared/sub-02/ses-/MRI/401-anat-T1w/header_dump_N2D_0001.json  
prepared/sub-02/ses-/MRI/401-anat-T1w/N2D_0001.nii
```

(Once per dataset)

- **prepared** – prepared dataset
- **bids** – place for bidsified dataset

bidsmap.yaml

1. *What is the purpose of this document?*

unknown.yaml

```

bids__ : 1.2.0
MRI:
  bidsmeNIFTI:
    __unknown__:
      - provenance: prepared/sub-02/ses-/MRI/4
01-anat-T1w/N2D_0001.nii
      checked: false
      suffix: ''
      attributes: {}
      bids: !!omap []
      json: !!omap []
      - provenance: prepared/sub-02/ses-/MRI/4
02-anat-T2w/N2D_0194.nii
      checked: false
      suffix: ''
      attributes: {}
      bids: !!omap []
      json: !!omap []

```

## Step 3b: Edit attributes

(Once per dataset)

bids/code/bidsme/bidsmap.yaml

bidsmap.yaml

```
bids__ : 1.2.0
MRI:
  bidsmeNIFTI:
    anat:
      - provenance: prepared/sub-02/ses-/MRI/40
1-anat-T1w/N2D_0001.nii
      checked: false
      suffix: 'T1w'
      attributes:
        ProtocolName: anat-T1w
        bids: !!omap []
        json: !!omap []
      - provenance: prepared/sub-02/ses-/MRI/40
2-anat-T2w/N2D_0194.nii
      checked: false
      suffix: 'T1w'
      attributes:
        ProtocolName: anat-T2w
        bids: !!omap []
        json: !!omap []
```

header\_dump\_N2D\_0001.json

```
"ManufacturerModelName": "nifti2dicom",
"PatientName": "Jane_Doe",
"PatientID": "02",
"PatientBirthDate": "1966-01-01",
"PatientSex": "F",
"OtherPatientIDs": "",
"PatientAge": 42,
"PatientWeight": 75.0,
"PregnancyStatus": 4,
"SliceThickness": 0.6666666666666667,
"SpacingBetweenSlices": 0.6666666666666667,
"SoftwareVersions": "0.4.11",
"ProtocolName": "anat-T1w",
"StudyInstanceUID": "1.2.826.0.1.3680043.2.1143.2592092611698916978113112155415165916",
"SeriesInstanceUID": "1.2.826.0.1.3680043.2.1143.515404396022363061013111326823367652",
"StudyID": "433724515",
"SeriesNumber": 401,
"AcquisitionNumber": 1,
```

# Step 3b: Edit bids entities

(Once per dataset)

## bidsme map prepared bids

### bidsmap.yaml

```
MRI:
  bidsmeNIFTI:
    anat:
      - provenance: prepared/sub-02/ses-/MRI/4
01-anat-T1w/N2D_0001.nii
      example: anat/sub-02_T1w
      checked: false
      suffix: 'T1w'
      attributes:
        ProtocolName: anat-T1w
      bids: !!omap
        - acq: <AcquisitionNumber>
        - ce: ~
        - rec: ~
        - run: <InstanceNumber>
        - mod: ~
      json: !!omap
        - ContrastBolusIngredient: ~
        - DeviceSerialNumber: '12345'
        - StationName: 'abcde'
        - MagneticFieldStrength: '7'
```

### header\_dump\_N2D\_0001.json

```
"PatientAge": 42,
"PatientWeight": 75.0,
"PregnancyStatus": 4,
"SliceThickness": 0.66666666534882,
"SpacingBetweenSlices": 0.66666666534882,
"SoftwareVersions": "0.4.11",
"ProtocolName": "anat-T1w",
"StudyInstanceUID": "1.2.826.0.1.3680043.2.1143.2592092611698916978113112155415165916",
"SeriesInstanceUID": "1.2.826.0.1.3680043.2.1143.515404396022363061013111326823367652",
"StudyID": "433724515",
"SeriesNumber": 401,
"AcquisitionNumber": 1,
"InstanceNumber": 1,
"PatientOrientation": [
  "L",
  "R"
],
"ImagePositionPatient": [
  -91.4495864331908,
```



## Step 3b: Check entries

(Once per dataset)

bidsme map prepared bids

bidsmap.yaml

```
bids__ : 1.2.0
MRI:
  bidsmeNIFTI:
    anat:
      - provenance: prepared/sub-02/ses-/MRI/401-anat-T1w/N2D_0001.nii
        example: anat/sub-02_acq-1_run-1_T1w
        checked: true
        suffix: 'T1w'
        attributes:
          ProtocolName: anat-T1w
        bids: !!omap
          - acq: <AcquisitionNumber>
          - ce: ~
```

- provenance – first matching file
- example – bidsified name
- checked – set to true if everything is ok

## Step 4: bidsify

### bidsme bidsify prepared/ bids/

```
bids/
|-- code/
|   |-- bidsme/
|   |   |-- .bidsmap.yaml.swp
|   |   |-- bidsify/
|   |   |   |-- bidsme.err
|   |   |   |-- bidsme.log
|   |   |-- bidsmap.yaml
|   |   |-- map/
|   |   |   |-- bidsme.err
|   |   |   |-- bidsme.log
|-- participants.json
|-- participants.tsv
|-- sub-02/
|   |-- anat/
|   |   |-- sub-02_acq-1_run-194_T1w.json
|   |   |-- sub-02_acq-1_run-194_T1w.nii
|   |   |-- sub-02_acq-1_run-1_T1w.json
|   |   |-- sub-02_acq-1_run-1_T1w.nii
|   |-- sub-02_scans.json
|   |-- sub-02_scans.tsv
```

File: sub-02\_acq-1\_run-1\_T1w.json

```
{
  "DeviceSerialNumber": "12345",
  "StationName": "abcde",
  "MagneticFieldStrength": "7",
  "Manufacturer": "BIOLAB",
  "ManufacturersModelName": "nifti2dicom",
  "SoftwareVersions": "0.4.11"
}
```

# Usefull links

- Code:
  - GitLab:  
`https://gitlab.uliege.be/CyclotronResearchCentre/Public/bidstools/bidsme/bidsme`
  - GitHub:  
`https://github.com/CyclotronResearchCentre/bidsme`
- Example dataset:
  - GitLab:  
`https://gitlab.uliege.be/CyclotronResearchCentre/Public/bidstools/bidsme/bidsme\_example`
  - GitHub:  
`https://github.com/CyclotronResearchCentre/bidsme\_examples`
- Contact:  
`Nikita.Beliy@uliege.be`