

QA/QC system for untargeted metabolomics: what we have learned from our previous studies


Anaïs RODRIGUES, Djulia BENSAADA, Pierre-Hugues STEFANUTO, Jean-François FOCANT

University of Liege (ULiege) – Chemistry, Sciences, MolSys Unit

ACS Fall 24 – August 18th, 2024 – Denver, USA



Introduction: metabolomics



Introduction: metabolomics

- Over the past decade, **metabolomics** - *the systematic study of the complete set of metabolites in a biological sample* - has gained increasing popularity and significance in the life sciences
- Alongside other “omics”, **metabolomics** provides **valuable insights** into **specific reactions** occurring in humans
- It helps us **understand various metabolic pathways** involved in **pathological processes**
- Metabolomics has been applied across **numerous fields**:
 - Environment
 - Toxicology
 - Disease diagnosis
 - *Etc.*



Introduction: QA/QC

- Quality assurance (QA) and quality control (QC) are two aspects of quality management
- QA: “part of quality management focused on providing confidence that quality requirements will be fulfilled”
- QC: “part of quality management focused on fulfilling quality requirements”
- Well known in industry and for targeted analysis
- **What about untargeted analysis and metabolomics?**



Quality System, Quality Assurance, and Quality Control Relationships

ISO 9000 – Quality management



Available online at www.sciencedirect.com



Journal of Chromatography A, 1019 (2003) 261–272

JOURNAL OF
CHROMATOGRAPHY A

www.elsevier.com/locate/chroma

Standardized test mixture for the characterization of comprehensive two-dimensional gas chromatography columns: the Phillips mix

Jean-Marie D. Dimandja^{a,*}, Garrick C. Clouden^b, Ivelisse Colón^c,
Jean-François Focant^d, Whitney V. Cabey^a, Ritchard C. Parry^e

^a Department of Chemistry, Spelman College, 350 Spelman Lane, SW Box 279 Atlanta, GA 30314, USA

^b Meharry Medical College, 1005 Dr. D. B. Todd Jr. Blvd., Nashville, TN 37208, USA

^c Pfizer Global Research and Development, Eastern Point Road, Groton, CT 06340, USA

^d Centers for Disease Control and Prevention, 4770 Buford Highway NE, Atlanta, GA 30041, USA

^e LECO Corporation, 3000 Lakeview Avenue, St. Joseph, MI 49085, USA

**analytical
chemistry**

pubs.acs.org/ac

Article

A Proposed Quality Control Standard Mixture and Its Uses for Evaluating Nontargeted and Suspect Screening LC/HR-MS Method Performance

Ann M. Knolhoff,^{*,†} Jacob H. Premo,[†] and Christine M. Fisher



Cite This: *Anal. Chem.* 2021, 93, 1596–1603



Read Online

Metabolomics (2022) 18: 70

<https://doi.org/10.1007/s11306-022-01926-3>

REVIEW ARTICLE




Quality assurance and quality control reporting in untargeted metabolic phenotyping: mQACC recommendations for analytical quality management

Jennifer A. Kirwan^{1,2,3} · Helen Gika^{4,5} · Richard D. Beger⁶ · Dan Bearden⁷ · Warwick B. Dunn⁸ · Royston Goodacre⁸ · Georgios Theodoridis^{9,5} · Michael Witting¹⁰ · Li-Rong Yu⁶ · Ian D. Wilson^{8,11} on behalf of the metabolomics Quality Assurance and Quality Control Consortium (mQACC)

GC×GC: instrumentation

LN2 Thermal Modulator



“Normal set” configuration


¹D: Rxi-5SilMS (30 m, 0.25mm ID, 0.25 μm)

²D: Rxi-17SilMS (2 m, 0.25mm ID, 0.25 μm)

Peltier Thermal Modulator




PEGASUS® BT 4D
GC×GC-TOFMS
LECO



GC oven Secondary oven Transfert line Modulator


Liquid injection (1μL)




Acquisition rate 200 Hz
Range: 35 – 450 mu



Compounds selection – QC37




Compounds selection – QC37




Cumulative QC charts – for instrument follow-up

- Building of a **cumulated** QC chart **within a year** with the **score** of all areas from the 37 compounds of our QC37




Cumulative QC charts – for instrument follow-up



- With PQN normalization and without internal standard (ISTD) correction

SD_{average}
45,8%



- With internal ISTD correction

SD_{average}
6,8%

QA for study follow-up

- Building of specific QC charts **for study follow-up**
 - Identification assurance, checked for carry over (Step 1)
 - Baseline for QC chart (Step 2)
 - Randomized sample analysis with QC samples (Step 3)
 - Repeating step 1 to check the system stability (Step 4)



International Journal of
Molecular Sciences




an Open Access Journal by MDPI

**Comprehensive Insight into Colorectal Cancer Metabolites and Lipids
for Human Serum: A Proof-of-Concept Study**

Kinjal Bhatt; Tiziana Orlando; Marie-Alice Meuwis; Edouard Louis; Pierre-Hugues Stefanuto;
Jean-François Focant

Int. J. Mol. Sci. 2023, Volume 24, Issue 11, 9614



QA for study follow-up



QA for study follow-up




MS Djulia
Bensaada






MS Djulia Bensaada

Baseline for QC Chart





Randomized bacterial sample analysis with QC samples



QA for study follow-up



MS Djulia Bensaada



Century mix project

- An interlaboratory worldwide project



Century mix project



- An interlaboratory worldwide project
- A mixture with ~100 compounds to characterize systems and column sets

Compounds	Function
Alkanes	Mark the direction of the primary column separation gradient. Hexane(C ₆) to Eicosane (C ₂₀) range.
Aromatic hydrocarbons	Mark the direction of the secondary column separation gradient. Benzene (1-ring) to Phenanthrene (3-ring) range.
Homologous compounds	Markers of the elution order in the secondary column separation gradient. Ten series of homologous compounds include ketones (C ₆ to C ₁₂), primary alcohols (C ₅ to C ₁₂), secondary alcohols (C ₅ to C ₁₀), FAMES (C ₅ to C ₁₃), alkyl benzenes (C ₇ to C ₁₅), alkenes (C ₈ to C ₁₁), alkyl acetates (C ₅ to C ₈), carboxylic acids (C ₄ to C ₈), Lactones (C ₆ to C ₉), and phthalates (C ₁₀ , C ₁₂ , and C ₁₄).
Grob mix and McReynolds compounds	Markers of selectivity metrics connected to 1D GC.
Miscellaneous compounds	Selected marker compounds connected to a variety of application areas (bioanalytical, environmental, foods and flavors)


Century mix project

- An interlaboratory worldwide project
- A mixture with ~100 compounds to characterize systems and column sets
- To be used as reference for QC in various applications




Century mix project

- Normal column set



Century mix project

- Reverse column set



Stay tuned for future work

- **Continue follow-up** on both Peltier and Cryo instruments and **implement QA/QC system** on our other instruments using both **QC37** and **CM**
- Development of **CM** as a standardized mixture expected to **enhance knowledge of column selectivity**
- **Further investigation of different column sets** planned to **gain deeper insight** into the relationship between **orthogonality and separation efficiency in GCxGC**
- **Create a retention index database** obtained under standard operating conditions across various stationary phases **to aid in column selection**



16TH Multidimensional Chromatography Workshop

 **LIÈGE, BELGIUM**

 **FEBRUARY 3-5, 2025**


FREE registration

PRESENTATIONS by key speakers

DISCUSSION on hot topics

POSTER SESSION with awards

www.multidimensionalchromatography.com



Liège University



Think tank



Networking hub

Acknowledgements



Prof. Jef Focant

Dr. PH Stefanuto

Dr. Kinjal Bhatt

Thibaut Dejong

Thibault Massenet

Anika Lokker

Géraldine Dumont

Djulia Bensaada



Dr. Armélinda Agnello

Stéphanie Hendrick

Grégory Gridelet

Vicky Borguet

Alex Haway

Océane Tirsell

Roberto ☐

Prof. John Dimandja
(FDA)

Dr. Giorgia Purcaro
(Gembloux Agro-Bio
Tech, University of Liège)

